



Ensure Patient Safety  
Protect your Investment

## DP850

Premium Steam Distillation Appliance

### The Best Steam-Distilled Water Solution for Dentists

**PROBLEM #1: BACTERIA & BIOFILM.** According to the American Dental Association, dental unit waterlines have been shown to harbor, in significant numbers, a wide variety of bacteria, fungi, and protozoans.

**PROBLEM #2: VERY SENSITIVE AUTOCLAVES.** New autoclaves and Statim machines require pure distilled water on a consistent basis or they will not operate properly and their lifespan could be shortened.

**SOLUTION:** Dentists should use only pure, steam-distilled water in their Dental Units and autoclaves.

### Recommended by Unthank Design Group

*"For patient safety, we have long advocated using steam-distilled water in a closed water system for dental delivery. A DENTAL PURITY in-office distillation system provides for a convenient continuous supply of the purest water."*



Michael Unthank DDS  
Registered Professional Architect & Dentist  
Unthank Design Group ([www.unthank.com](http://www.unthank.com))

## DP850

### Premium Steam Distillation Appliance

The DP850 Steam-Distillation Appliance is an extremely reliable and convenient appliance for dental offices. The DP850 automatically produces up to 7 gallons of steam-distilled water each day and has a 5 gallon storage tank. Put this appliance on a shelf, or purchase the optional stand and simply dispense the water from the front of the appliance when you need it.

- Automatic operation.
- Produces plenty of steam-distilled water for small to large dental offices.
- Simple installation.



Visit us on-line at [www.DentalPurity.com](http://www.DentalPurity.com)

Or call one of our customer care representatives today at **1-800-617-2401.**

## PROVIDE SAFE, PURE WATER FOR YOUR PATIENTS

Steam-distilled water is an important key in preventing bacteria and biofilm build-up in your water lines so you can provide a safe, pure water supply for your patients.

## TAKE CONTROL OF YOUR DISTILLED WATER

The DP850 allows you to take control of your own supply of pure, steam-distilled water. No more running out of distilled water, last minute trips to the store and no more of the hassle or wasted space of bottled water.

## AUTOMATIC OPERATION

The DP850 automatically produces plenty of steam-distilled water each day. Simply dispense the distilled water directly from the front of the DP850 whenever it's needed.

## 40 YEARS OF TRUST

For 40 years, dentists around the world have used our premium steam distillation appliances (under different brand names) to ensure that they are using the highest quality water at all times.

## FAR SUPERIOR TO COMPETITIVE DISTILLERS

The DP850 is far superior to competitive distillers on the market today. Here's why...

**QUALITY CONSTRUCTION** The DP850 is constructed of stainless steel and top-quality components to provide years of dependable service.

**AN EXCELLENT INVESTMENT** When comparing the DP850 to competitive distillers, which have a very short lifespan, or bottled water it's clear that the DP850 is an excellent investment that will pay you back many times over.

**AUTOMATIC OPERATION** The DP850 automatically produces distilled water and stores it in a 5 gallon internal stainless steel storage tank so you'll never run out of distilled water. Just dispense distilled water whenever it's needed.

**FAR SUPERIOR TO PLASTIC DISTILLERS** The DP850 is superior to competitive distillers in design and construction. It's automatic operation offers more convenience and superior water quality.



## SPECIFICATIONS

item #.....	DP850
width.....	21" (53.3 cm)
depth.....	18" (45.7 cm)
height.....	16.5" (41.9 cm)
height (w/floor stand).....	34.5" (87.6 cm)
dry weight.....	48 lbs. (18.9 kg)
current draw (120V).....	9 amps
current draw (240V).....	4.5 amps
total wattage.....	1,100 watts
daily production.....	7 gallons
storage capacity.....	approx. 5 gal.

Specifications may be subject to change without notice.

**CONSISTENT PURITY** The DP850 provides unsurpassed purity on a consistent basis.

**LONG LIFESPAN** The DP850 is a heavy duty appliance and is built-to-last. The DP850 has a proven track record and a 15-year limited warranty.

**FREE DRINKING WATER** Save even more money each month by filling up bottles of drinking water with fresh, great-tasting steam-distilled water.

**SIMPLE INSTALLATION** The DP850 is very simple to install. It just requires a 1/4" line into the machine. No drain line is required.

**EASY MAINTENANCE** The DP850 requires very little maintenance.