



*Ensure Patient Safety
Protect your Investment*

Have you ever had this experience...

You're not getting water out of your high-speed hand piece so you call a technician and they find green or brown slime in the tubing. This is biofilm. Biofilm is different types of bacteria living together in a colony that is strong and hard to eliminate.

Biofilm is an important patient safety issue that the ADA wants all dentists to take seriously. Prevention, by only using pure, distilled water is vital.

Or maybe you've had this experience...

Your autoclave needs frequent cleaning or simply won't function because the water that you are using isn't as pure as it needs to be. Or perhaps you've actually ruined your autoclave because of poor quality water. Again, the solution is to use only distilled water!

All dental offices need a reliable source of distilled water in order to ensure the highest level of patient safety and to protect their expensive dental equipment.

Visit us at www.DentalPurity.com

Or call us at **1-800-617-2401**

prevent & protect

biofilm build-up in your dental units

Biofilm build-up in your dental units is a serious patient safety issue that can be a legal risk and can affect patient confidence. The ADA encourages all dentists to take this issue seriously. By only using distilled water in your dental units you are taking an important step to prevent bacterial and biofilm contamination.

your expensive dental equipment

Common contaminants in tap water such as hardness, minerals, and algae can interfere with your expensive dental equipment, shorten the lifespan and greatly increase cleaning and maintenance frequency. By using only pure, distilled water in your dental equipment, such as autoclaves, you can save thousands of dollars over the years.

Distilled water is water that has been boiled, sterilized and purified. All dental offices need a reliable source of distilled water in order to ensure the highest level of patient safety and to protect their expensive dental equipment.

Perhaps you're already using distilled water. That's good, but if you are using bottled water or a low-quality distiller you are most likely spending way too much money. After all, bottled water costs more than gasoline! Dental Purity brand appliances offer dentists the most cost effective solution for distilled water.

DP360

\$995 15-Year Cost
15 Years Plus Lifespan



Value Comparison

Competitive Distillers

\$3,600 15-Year Cost
1-2 Year Lifespan



Bottled Water

\$6,300 15-Year Cost



FORMULA: $\$995 \times 1$ (15 year lifespan) = $\$995$; $\$480 \times 7.5$ (2 year lifespan) = $\$3,600$; $\$35$ per month $\times 15$ years = $\$6,300$.